


Roll-a-Ramp®

Where strength and flexibility meet!



Roll-a-Ramp®

Strength, flexibility and portability come together in a Roll-A-Ramp®. They can easily be lengthened, shortened and rolled up for transport and storage. Roll-A-Ramp® Modular Ramp Systems can be used anywhere access is limited.

The Roll-A-Ramp® system gives you the versatility of still having a portable system — at any length! Roll-A-Ramp® is the most portable, versatile ramp available. It is great for wheelchair and scooter users who want an alternative to chair lifts, permanent ramps or van conversions. It is also great for anyone looking for a superior quality ramp for any use.

FEATURES INCLUDE:

- Lightweight and Strong — 6061 rated Aircraft Aluminium
- Portable — Take it with you
- Versatile — Any size, any time.
- Unique Design — Patented link construction means you will never have the wrong size ramp; simply add additional links to change length with simple tools provided with your ramp. If your situation changes? Change the ramp to meet your needs.

WIDTH AVAILABLE IN

300mm (twin track), 760mm and 915mm wide

LENGTH

These ramps can be supplied to any length you require. The following are two common size examples:

Example technical specifications for 760mm wide x 1800mm long ramp

Approx weight: 14.0kg
Safe working load: 454kg

Example technical specifications for 915mm wide x 1800mm long ramp

Approx weight: 16.33kg
Safe working load: 454kg

Custom sizes available - let us know the length you need. Please enquire.

300mm wide ramp



Optional removable handrails

Available as an addition to the Roll-A-Ramp. All aluminium construction. Heavy duty and lightweight. Protective coating. Never rust or corrode. Easily installed or removed.



Optional removable handrails

Available as an addition to the Roll-A-Ramp. All aluminium construction. Heavy duty and lightweight. Protective coating. Never rust or corrode. Easily installed or removed.



Optional support stands

The support stands are recommended for ramps over 3700mm in length. Ideal for situations where the ground surface is uneven or extra support is required. Each leg may be independently adjusted.



Additional Links

Use additional links to lengthen the ramp. Links easily attach with simple tools. Never have the wrong size ramp.



Quick Release Pins

Easily separate your Roll-A-Ramp into shorter, even lighter sections for easy storage, handling and transport.



Roll-A-Ramps can be used in almost any situation

CUBRO

healthcare equipment solutions

Head Office: 149 Taurikura Drive, Tauriko, Tauranga 3171 | PO Box 9144, Greerton, Tauranga 3142, New Zealand
Phone: +64 7 578 722 8 | Freephone: 0800 656 527 | Email: sales@cubro.co.nz | cubro.co.nz